

# Selectin Biosciences Inc.

## $\alpha$ (2-3, 6, 8, 9) Neuraminidase

$\alpha$  (2-3, 6, 8, 9) neuraminidase (EC 3.2.1.18) cleaves all non-reducing terminal sialic acid residues from complex carbohydrates and glycoproteins. The relative cleavage rates for different linkages are:  $\alpha$ (2-6) >  $\alpha$ (2-3) >  $\alpha$ (2-8),  $\alpha$ (2-9).

In addition, the enzyme will cleave branched sialic acids (linked to an internal residue). This property makes it unique among sialidases (Figure 1). High concentrations of enzymes and prolonged incubation times may be required for cleaving branched residues. To cleave only non-reducing terminal  $\alpha$ (2-3) unbranched sialic acid residues, use  $\alpha$  (2-3) neuraminidase(GE20).  $\alpha$  (2-3, 6, 8, 9) neuraminidase is isolated from a clone of *Arthrobacter ureafaciens*. The enzyme has been extensively characterized using oligosaccharide standards.

$\alpha$  (2-3, 6, 8, 9) neuraminidase is useful for:

- Structural analysis of oligosaccharides
- Determining sialic acid linkage
- Glycoprotein deglycosylation
- Removing heterogeneity from glycoproteins

**Product Code:** GE 23

## Specifications

**Activity:**  $\geq$ 135 U/mg, 5 U/mL

**Storage:** Store at 4°C. Do not freeze.

**Formulation:** The enzyme is provided as a sterile solution in 20 mM Tris pH 7.5, 25 mM sodium chloride.

**Stability:** Stable at least 12 months when stored properly. Several days exposure to ambient temperatures will not reduce activity.

## Product Description

**Molecular weight:** 70,000 Daltons

**Purity:** Each lot of  $\alpha$ (2-3, 6, 8, 9) neuraminidase is tested for contaminating protease as follows; 10  $\mu$ g of denatured BSA is incubated for 24 hours with 2  $\mu$ L of enzyme. SDS-PAGE analysis of the treated BSA shows no evidence of degradation.

The production host strain has been extensively tested and does not produce any detectable glycosidases.

**Specificity:** All non-reducing terminal branched and unbranched sialic acids (see Figure 1).

**pH Range:** Optimum: pH 6  
Range: pH 4.5 - 7

The recommended buffer concentrate provides the optimal pH for enzyme activity with the standard substrate. If glycosidase treatment is performed at suboptimal pH because of glycoprotein solubility or activity requirements, expect some diminution in enzyme activity.

## Assay

One unit of  $\alpha$  (2-3, 6, 8, 9) neuraminidase activity is defined as the amount of enzyme required to produce 1  $\mu$ mole of methylumbelliferone in 1 minute at 37°C, pH 5

# Selectin Biosciences Inc.

from 2'-(4-methylumbelliferyl)- $\alpha$ -D-N-acetylneuraminic acid.

## Reagents

- 5X Reaction buffer 6 – 250 mM sodium phosphate pH 6.0

## Suggestions For Use

### Procedure for Desialylation

Desialylation may be monitored by SDS-PAGE if the size differential between native and desialylated protein is sufficient for detection.

1. Add up to 100  $\mu$ g of glycoprotein or 1 nmol of oligosaccharide to tube.
2. Add water to 13  $\mu$ L and 4  $\mu$ L 5X Reaction Buffer.
3. Add 2  $\mu$ L  $\alpha$  (2-3, 6, 8, 9) neuraminidase.
4. Incubate at 37°C for 1 hour.

NOTE: longer incubation times are necessary if branched sialic acids are present.

## References

Cornfield, A. P., H. Higa, J. C. Paulson and R. Schauer. The specificity of viral and bacterial sialidases for alpha(2-3) and alpha(2-6)-linked sialic acids in glycoproteins. **Biochim Biophys Acta** **744**:121-126 (1983).

Dwek, R. A. , C. J. Edge, D. J. Harvey, M. R. Wormald and R. B. Parekh. Analysis of glycoprotein-associated oligosaccharides. **Ann Rev Biochem** **62**:65-100 (1993).

Kobata, A. Use of endo- and exoglycosidases for structural studies of glycoconjugates. **Anal Biochem** **100**:1-14 (1979).

Ohta, Y., Y. Tsukada and T. Sugimori. Purification and properties of neuraminidase isoenzymes in *Arthrobacter ureafaciens* mutant. **J Biochem (Tokyo)** **106**:1086-1089 (1989).

Prime, S. J. Dearnley, A. M. Venton, R. B. Parekh and C. J. Edge. Oligosaccharide sequencing based on exo- and endoglycosidase digestion and liquid chromatographic analysis of the products. **J Chromatogr A** **720**:263-274 (1996).

Uchida, Y., Y. Tsukada and T. Sugimori. Enzymatic properties of neuraminidases from *Arthrobacter ureafaciens*. **J Biochem (Tokyo)** **86**:573-585 (1979).

Selectin Biosciences Inc warrants that the above product conforms to the specifications described herein. Should the product fail for reasons other than through misuse Selectin Biosciences Inc. will, at its option, replace free of charge or refund the purchase price. This warranty is exclusive and Selectin Biosciences Inc. makes no other warrants, expressed or implied, including any implied conditions or warranties of merchantability or fitness for any particular purpose. Selectin Biosciences Inc. shall not be liable for any incidental, consequential or contingent damages.


This product is intended for in vitro research only.

REVISION 5/22/12

# Selectin Biosciences Inc.

**Figure 1**

**Linkage specificities showing cleavable residues (in bold) for  $\alpha(2-3, 6, 8, 9)$  Neuraminidase**



Gal = Galactose; Glc = Glucose; Man = Mannose; GalNAc = N-acetylgalactosamine; GlcNAc = N-acetylglucosamine; NeuAc = N-acetylneuraminic Acid (Sialic Acid)